

SPECTRA

CARBON-BRUSH MICROMOTOR

B-350

DIRECT CONTROL. EVERY DETAIL.



35,000

MAX. RPM

CARBON

MOTOR TECHNOLOGY

ON / OFF

FOOT CONTROL

DIGITAL SPEED DISPLAY

Clear speed indication up to 35,000 RPM.

HAND / FOOT SELECTION

Set speed at the dial; start/stop with the pedal.

FORWARD / REVERSE

Choose rotation direction for the task.

2.35 mm SHANK

Designed for the specified matching bur shank.

SPECTRA CERAM TRADING L.L.C.

Office 505, City Bay Business Centre, Hor Al Anz East, Dubai, UAE

www.spectra-ceram.com
info@spectra-ceram.com
+971 56 979 7090