

SPECTRA

BRUSHLESS MICROMOTOR

BL-500

MORE RANGE. VARIABLE CONTROL.



50,000

MAX. RPM

BRUSHLESS

MOTOR TECHNOLOGY

VARIABLE

FOOT CONTROL

BRUSHLESS DRIVE

No carbon-brush replacement in the motor.

VARIABLE FOOT CONTROL

Pedal pressure adjusts speed within the dial limit.

DIGITAL & REVERSIBLE

0-50,000 RPM display range; FWD / REV selection.

SETTING MEMORY

Remembers speed, direction and control mode.

SPECTRA CERAM TRADING L.L.C.

Office 505, City Bay Business Centre, Hor Al Anz East, Dubai, UAE

www.spectra-ceram.com
info@spectra-ceram.com
+971 56 979 7090